

1ST-COMFORT™

TROLLEY MATTRESS



- HIGH COMFORT MATTRESS
- OPTIMIZED FOR X-RAY
- HYGIENIC AND DURABLE



1st-Comfort™ is a pressure-relieving mattress. Esongrip™ is easy to adjust and comes with adjustable head cushion.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUKTINFORMATION

The mattress contains a 20mm pressure-relieving layer in viscoelastic foam, volume-weight 52kg/m³. The core is made of 60 mm:s flexible foam, volume-weight 38kg/m³. Both materials is compatible to the cigarett- and ignition test according to Swedish standard, without any addition of flame retardents. The fabric is resistant to urin exposure and cleaning with alcohol. Production without use of CFC. 1st-Comfort comes with a 1 year technical guarantee of function.

ART.NO	NAME
1STCOM	1st-Comfort trolley mattress Black (complete with cushion) 60 x 195 x 8cm
1STBAS	1st-Basic trolley mattress Black (without grip and low friction underside) 66 x 203 x 7cm
1STKUDDE	1st-Comfort head cushion Blue 50 x 22 x 8cm

Accessories/spare parts

Special shape or size.

Washing

The fabric is cleaned with soap and alcohol.

CE Registrationcode 344-2757-35330

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"