

BEDSLIDE™ BS

TRANSFERRING / LIFTING



- ERGONOMIC
- HYGIENIC AND DURABLE
- EASY TO APPLY

BedSlide™BS is a flexible sliding board used in the transfer and turning of patients. It has been proven in use in operation-rooms, x-ray departments, emergency rooms, intensive care units and normal care units.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUKTINFORMATION

BedSlide™BS is made of 3 mm HD-Polyethene. The material is very flexible and smooth. Two handles allows for an ergonomic grip. BedSlide™BS contains it's flexibility even in tough pre-hospital conditions during winter. BedSlide™BS gives no artefacts in X-ray. BedSlide™BS comes with a warranty regarding function 1 year.

ART.NO
BEDSLIDE
BEDSLIDEBS
BEDSLIDEL
BEDSLIDEXL
BEDSLIDASET

NAME
BedSlide sliding board 37 x 80cm
BedSlideBS sliding board 37 x 80cm
BedSlide Large sliding board 37 x 100cm
BedSlide Extra Large sliding board 37 x 80cm
8pc. BedSlide, 4pc. BedSlide L och 2pc. BedSlide XL

Extras/spare parts

Special shape, special size, safety belts.

Care

The material does not absorb any fluid and is easy to rinse. Cleaning with soap-alcohol solution.

CE Registrationcode 344-2757-35328

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"