


# BEDSLIDE™

## TRANSFERRING / LIFTING



- ERGONOMIC
- HYGIENIC AND DURABLE
- EASY TO APPLY

BedSlide™ is a flexible sliding board used in the transfer and turning of patients. It has been proven in use in operation-rooms, x-ray departments, emergency rooms, intensive care units and normal care units.



ESON COMFORT  
Idé-Möbler i Landeryd AB  
Centrumvägen 7  
SE-314 94 LANDERYD

Tel +46 (0)371-40051  
[info@esoncomfort.com](mailto:info@esoncomfort.com)

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

### PRODUKTINFORMATION

BedSlide™ is made of 3 mm HD-Polyethene. The material is very flexible and smooth. Two large handles allows for an ergonomic grip. BedSlide™ contains it's flexibility even in tough pre-hospital conditions during winter. BedSlide™ gives no artefacts in X-ray. BedSlide™ comes with a warranty regarding function 1 year.

ART.NO	NAME
BEDSLIDE	BedSlide sliding board 37 x 80cm
BEDSLIDEBS	BedSlideBS sliding board 37 x 80cm
BEDSLIDEL	BedSlide Large sliding board 37 x 100cm
BEDSLIDEXL	BedSlide Extra Large sliding board 37 x 120cm
BEDSLIDASET	8pc. BedSlide, 4pc. BedSlide L och 2pc. BedSlide XL

#### Extras/spare parts

Special shape, special size, safety belts.

#### Care

The material does not absorb any fluid and is easy to rinse. Cleaning with soap-alcohol solution.

CE Registrationcode 344-2757-35328

ESON COMFORT  
www.esoncomfort.com  
"WE INVENT THE FUTURE"