

NEFROLIFT™

POSITIONING



- COMFORTABLE
- ERGONOMIC
- HYGIENIC

Nefrolift™ is a lifting device and aid used in patient positioning during kidney diagnostics, ultra sound guided kidney punction or biopsies.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

Nefrolift™ has a padding consisting of a 15 mm pressure relieving polyurethane foam with a volume weight of about 52kg / m³. Inflatable inner part with the dimensions 40x30cm and a hand pump with valve. Under the padding and the bladder is an acrylic glass plate to stabilize the lift movement. Nefrolift™ has dimensions 50x40x2,5cm at rest and 50x30x20cm when activated. The cover is made of syntetic leather and can handle cigarette and ignition tests in accordance with Swedish minimum standard without the addition of flame retardants. Production takes place without CFC. Nefrolift™ has a technical function warranty of 1 year.

ART.NO	NAME
NEFRO	Nefrolift 40 x 50 x 2,5 cm, with hand pump

Special shape or size. Accessories/spare parts

Washing
The cover is cleaned with soap and alcohol.

CE Registrationcode 344-2757-81424

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"