


TT-COVER K™ MINI

TRAUMA



- WARMING
- HYGIENIC
- TIMESAVING



TT-Cover K™ Mini is a disposable cover with a separate soft/hard head cushion. It is designed for optimum tensile strength. TT-Cover K™ Mini is for paediatric use, 1 - 8 years.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Cover K™ Mini is a disposable cover with a separate soft/hard head cushion. It is developed to achieve durability and safe immobilisation. The fabrics is composed to reduce the loss of patientwarmth and to prevent contamination of the TraumaTransfer™-mattress and the trolley. It is designed for optimum tensile strength. Tested for a work load of 150 kg in wet condition. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic wast when burning.

ART.NO	NAME
TTCOVERK	TT-Cover™ disposable cover for TraumaTransfer™ (25 pcs)
TTCOVERKM	TT-Cover™Mini disposable cover for TraumaTransfer™ Mini (25 pcs)
TTHEADFIXH	Head cushion for TT-Cover K™ (25 pcs) (order separately)
TTHEADFIXMH	Head cushion for TT-Cover K™ Mini (25 pcs) (order separately)

Extras / Spareparts

Special size, special shape.

CE Registrationcode 344-2757-118173

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"