



SPACE™

POSITIONING



- ERGONOMIC
- OPTIMIZED FOR X-RAY
- HYGIENIC AND DURABLE



SPACE™ is the leading system for patient positioning. The image above shows a complete set of SPACE™. On the right is a special wedge and SPF13 - Finger Stairs.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

Solid EVA with closed cells with a volume weight of 50kg/m³. High friction material. Gives no artefacts. No CFC. Space™ comes with a technical guarantee of function 1 year.

ART.NO	NAME	ART.NO	NAME
SPL38G	Space Large Gavel 38 cm	SPL50G	Space Large Gavel 50 cm
SPL60G	Space Large Gavel 60 cm	SPF38G	Space Flack Gavel 38 cm
SPS30	Space Small Small Knee 30 cm	SPM18R	Space Mini radius 18cm
SPF13	Space Hand profile Finger 13 cm	SPKS19	Space Wedge Small 0,5-4 cm
SPKL19	Space Wedge Large 0,5-6 cm	SP2040	Plate for X-ray 20 x 40 x 5 cm
SP1	SET av cushions (the large image)		

Accessories/spare parts

Special shape or size.

Washing

The fabric is cleaned with soap and alcohol.

CE Registrationcode 344-2757-35329

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"