

TT-COVER 2014™ MINI

TRAUMA



- WARMING
- HYGIENIC
- TIMESAVING



TT-Cover™ Mini is a disposable cover with an integrated head support. It is designed for optimum tensile strength. The 2014 version is upgraded in many ways and has been made a lot more user friendly. A thinner fabric has been developed with the same strength as previous version. For paediatric use, 1 - 8 year.



ESON COMFORT
 Idé-Möbler i Landeryd AB
 Centrumvägen 7
 SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Cover™ Mini is a disposable cover with an integrated head support. It is developed to achieve durability and safe immobilisation. The fabric is composed to reduce the loss of patient warmth and to prevent contamination of the TraumaTransfer™ Mini-mattress and the trolley. It is designed for optimum tensile strength. Tested for a work load of 150 kg in wet condition. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic waste when burning.

| ART.NO | NAME |
|-----------|---|
| TTCOVER2 | TT-Cover™ disposable cover for TraumaTransfer™ (25 pcs) |
| TTCOVER2M | TT-Cover™ Mini disposable cover for TraumaTransfer™ Mini (25 pcs) |

Extras / Spareparts

Special size, special shape.

Swedish patent search: 0402834-6

Publication: 529895IPC-klass: A61G 1/04, A61F 5/37

CE Registrationcode 344-2757-75049UNSPSC Code: 42.17.17.01GMDN Code: 35181

ESON COMFORT
www.esoncomfort.com
 "WE INVENT THE FUTURE"